# A-Series

Under vehicle lift

SV	5	Installationsanvisning / Bruksanvisning
EN	25	Installation instructions / User manual
DE	45	Einbauanleitung / Benutzerhandbuch
FR	65	Manuel d'installation / Manuel d'utilisation
ES	85	Instrucciones de instalación / Modo de empleo
IT	105	Istruzioni per l'installazione / Manuale delle istruzioni

SV

EN

DE

FR

ES

IT

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Thank you for choosing a

# A-Series from BraunAbility!

The manual that follows is an integral, important part of the product, providing you with information on how to perform accurate and safe installation of the product. Please read through the manual carefully before starting your work.

If you have any questions about your customer's adaptation, please feel free to contact us. Your feedback is important to us. It helps us in our work of continuously improving our product and its installation. You can find more information about us and our range of car adaptation solutions by visiting [www.braunability.eu](http://www.braunability.eu).

Once again, thank you for placing your confidence in our products!

## Installation

---

26	Pallet contents on delivery
27	Safety
28	Preparing the cassette for installation
29	Installing the cassette
30	Adjusting the outer end position
31	Install the pump module
32	Routing hose and cables
34	Routing and connecting the electrical supply
36	Adjusting the upper limit
37	Adjusting the inclination
38	Verification by the installer

## Service check list

---

40

## Technical specification


---

43


Diagrams, exploded views and manufacturer certificate (Product certificate) at the end of the manual.

# Installation



## Pallet contents on delivery

x 1 Lift cassette


x 1 Pump module with cover

x 1 Pump module handle

x 1 Hand-held control

x 1 Holder

x 2 Clips

x 5 Screws

x 1 Long red cable


x 1 Short red cable


x 1 Black cable

x 1 Extra cable lug

x 1 Plastic protector

x 1 Fuse holder

x 1 Fuse

x 2 Screws

## Safety

This lift has only been tested and approved for transferring people with reduced mobility into a vehicle.


The lift must not be used as a link-bridge to the pavement or raised groundwork. In order to achieve optimum embarkation and alighting, the lift platform should be flat on the ground.

Wheelchair brakes should always be applied during the lifting operation. The front roll stop and the inner barrier, see illustration below, should be raised during operation and be clear of any interference throughout. This is especially important to consider when an attendant is on the platform.

A safety valve is triggered when the maximum weight capacity of 350 kg is exceeded, when this is activated the platform will not lift. Should this happen, reduce the load and try again.


## Preparing the cassette for installation


*Always use customary personal protective equipment during installation.*

- 1 Raise the vehicle to facilitate access during the entire installation. Start by mounting the brackets, with the accompanying fasteners, delivered in a separate package. Follow the instructions supplied with the brackets. Note that the brackets should only be tightened lightly so there is some free play for adjusting the positioning of the cassette.
- 2 Cut away the banding around the cardboard packaging. Remove the cardboard hood, the boxes and the wooden transport brackets. Leave the cassette on the pallet.
- 3 Temporarily remove the screws and associated washers from the four mounting slides located along the sides of the cassette.

## Installing the cassette


- 1 Raise the pallet using a forklift. Position the cassette between the brackets, at least 10 mm from the rear axle.
- 2 Adjust the mounting slides so they correspond with the respective brackets. Fasten the cassette to the brackets using the screws and washers removed in an earlier step. Use the rearmost two holes. When attaching a footstep, an accessory that can be ordered separately, use the same mounting slide. The third hole is used only in this instance.
- 3 Lower the pallet and remove it. Retighten the screws to 50 Nm. Then retighten all the bracket fastenings onto the vehicle chassis.

## Adjusting the outer end position


If the Cassette lift is installed in combination with a footstep the outer end position of the lift needs to be adjusted. This to ensure that the bridge plate reaches sufficiently far into the vehicle. Adjust the outer end position before fitting the footrest.

- 1 Fold down the front cover and fix it in place with the wing screw.
- 2 Move the metal tab closer to the switch using the secondary hole pattern. Place the two screws in center position of the slots.

## Install the pump module


- 1 Install the pump module in the desired location inside the vehicle. Remember to place the module on the nearside of the vehicle to facilitate easier manual operation in the event of a lift failure. Use a suitable fastening in each of the three mounting points on the floor and in at least one of the mounting points on the side panel. Ensure the installation is sufficiently stable to withstand operation in an emergency when using the manual hand pump.
- 2 Store the handle mentioned above close to the pump module so it is easily accessible when required. Use the accompanying clips.
- 3 Remove the yellow transport plug from the oil reservoir air hose and replace it with the supplied vent cap.

**TIP!**  
Using the pump module cover as a template, cut out a corresponding opening in the vehicle trim panel. In this way, the pump module can be attached directly to the chassis and will protrude only a few centimetres from the panel.

## Routing hose and cables

- 1 Hydraulic hoses and electrical couplings between the power pack and the cassette, can be routed through the same hole. Use existing holes and spaces in the vehicle chassis for this purpose if possible. If this is not possible, drill a hole in a suitable position with a 50 mm hole saw. Place the plastic grommet into the hole to protect the Cables against sharp edges.
- 2 Lay the oil hose from the pump module to the cassette in an appropriate manner. Note that bends in the hose must not be less than 30 mm in diameter in order to ensure correct function. Connect the hose to the quick-coupling on the cassette.


- 3 Lay the wiring loom from the cassette to the pump module in an appropriate manner. The connectors are marked A and B. Plug into the corresponding sockets on the pump module.
- 4 When the lift is fully inserted into the cassette, wire C36 gives a + (plus) 12V signal. If a - (minus) 12V signal is desired via the vehicle ignition, see the connection diagram at the end of this manual.

Note! Cable C36 is only intended for connecting the lamp.

- 5 Connector to the Interlock kit (option). See separate instruction with the kit.

## Routing and connecting the electrical supply


- 1 Fit the fuse holder along with fuse in close proximity to the vehicle battery.
- 2 Connect the long red cable with the angled end lug to solenoid in the pump module. Route the cable towards the battery through the drilled hole in the vehicle floor and connect it to the fuse holder and make connection secure. The length of the long red cable can be shortened if too long. In this case, use the supplied cable lug that fits on the fuse holder.


Connect the short red cable. Mount one end to the fuse and the other end up to the 12V positive terminal on the vehicle battery.

- 3 Connect the black ground cable. Fasten one end to the pump module and the other end to a suitable grounding point. A suitable grounding point is an area of clean metal on the vehicle chassis, free from dirt and paint.
- 4 Connect the hand-held control to the corresponding socket on the pump module. Fasten the cable with the p clip to minimise strain on the connection. Fit the holder for the hand-held control in close proximity to the lift so that it is easily accessible when needed.

See supplementary illustration on next page.

### Pump module


When electric supply is connected a complete test run of the lift's full movement cycle must be performed. The lift will stop automatically once it has reached its set limits.

Keep a close eye on the movements of the lift during the cycle test and make adjustments according to the following pages if needed. Neglecting to do this could result in physical damage on vehicle and/or lift.

## Adjusting the upper limit

- The inner platform (bridge plate) should be level with the floor. If necessary, adjust the upper limit as illustrated below.

### Lift stops too high


Tap the clamp ring lightly to move it closer to the microswitch. The closer the clamp ring is to the microswitch, the lower the lift will stop. Tighten the screws to fix the new position.

### Lift stops too low

Tap the clamp ring lightly to move it further away from the microswitch. The further the clamp ring is from the microswitch, the higher the lift will stop. Tighten the screws to fix the new position.

## Adjusting the inclination

- 1 The inner platform (bridge plate) should be level with the floor. If necessary, adjust the inclination as illustrated below.
- 2 Adjust the cam. Ensure the wheel is straight in relation to the cam. Failure to do so can cause damage the bridge plate.
- 3 Manually deploy the travelling frame, see separate emergency instruction, in order to ensure the platform moves freely in an out of the cassette. If it does not, it may be because the mounting brackets are over tightened and may have deformed the sides of the cassette. The brackets should be adjusted if this is the case.

## Verification by the installer

According the directive EN 1756-2:2004+A1:2009, the installer shall confirm the compatibility between the lift and the vehicle.

Static test.....

### Deformation

- The unladen platform is positioned mid-way between ground level and vehicle floor level and measurements are taken of the height of the platform and its angular attitude relative to the vehicle floor.
- An evenly spread load of 125 % (438 kg) of the maximum load (350 kg) is applied to the platform, and subsequently removed.
- By repeating measurements of the height and attitude of the platform, verify that no permanent deformation has occurred in any part of the tail lift or its attachment to the vehicle which could affect the function of the tail lift.

### Drift

- An evenly spread load of 125 % (438 kg) of the maximum load is applied to the platform, positioned at floor vehicle floor level!. Measurements are taken of the height of the platform and its angular attitude relative to the vehicle floor. These measurements are repeated after a 15 min test period.
- Verify that the vertical drift of the platform between the two measurements has not exceeded 15 mm.
- Verify that the angular drift of the platform between the two measurements has not exceeded 2°.

Dynamic test.....

- With the maximum load applied to the platform, verify that the tail lift is able to operate throughout its full range of normal lifting, lowering and any tilting movements. If the lift is unable to lift with this weight, the lifting capacity must be increased. Proceed as follows (see also the illustration on the next page):
  1. Loosen the locking nut. Use a 14 mm (or 9/16") socket or swivel joint wrench.
  2. Turn the adjusting screw in a clockwise direction. Use a 3/16" Allen key.
  3. Retighten the locking nut. Use the Allen key to hold the adjusting screw in position when tightening the locking nut.

Test of operation and safety functions.....

- All functions of the tail lift and operations of all safety devices are verified after the static and dynamic tests have been completed. These tests do not apply to pipe break valves nor non-resettable safety devices such as electrical fuses (These items are the subject of a manufacturer's type test).

Test to verify that the tail lift cannot.....  
lift excessive load

- A load of 438 kg (125 % of the maximum load) is applied to the platform, positioned at ground level. Actuate the "UP" control and verify that the platform does not lift (tilt is permissible). If the lift is able to lift with this weight, the lifting capacity must be decreased. Proceed as follows (see also the illustration below):

1. Undo the locking nut. Use a 14 mm (or 9/16") socket or swivel joint wrench.
2. Turn the adjusting screw in a counterclockwise direction. Use a 3/16" Allen key.
3. Retighten the locking nut. Use the Allen key to hold the adjusting screw in position when tightening the locking nut.


Lift model \_\_\_\_\_

Serial number \_\_\_\_\_

Perform date \_\_\_\_\_

\_\_\_\_\_  
Name of dealer

\_\_\_\_\_  
Signatur, installation engineer

# Service check list

To ensure maximum service life and functionality please follow the check list below.

Lift model _____	Dealer _____
Serial number _____	Resp. mechanic _____
Vehicle reg. number _____	Owner _____
Vehicle make/model _____	Owner address _____




\*Hydraulic oil, BraunAbility art. no. 440612

\*\*OKS 375I, BraunAbility art. no. 400878



\*\*\*Light oil








EVERY DAY	A	Check that the bridge plate (A) goes up and down in the manner intended.
	B	Check that the hand rails with protective plates (B) are functioning in the manner intended. Protective plates are an option.
	C	Check that the outer roll stop (C) goes up and down in the manner intended.
	D	Check that the cassette (D) is free of external damage.
	E	Check that the handle (E) for the pump module is in place.
	F	Check that the hand control (F) and cable are intact and functioning.
		Check that there is no visible leaking of hydraulic oil.
		Check that instructions are available and in a legible condition.
EVERY MONTH	G	Check that the level of the *hydraulic oil is correct. When the lift is fully inside the cassette, the level must lie between the "stow level" mark and up to a maximum of 5 mm below the mark (G).
	H	Check that the fuses (H) are intact.
		Check that visible cables and tubing are free of external damage.
		Check that everything works in the manner intended.
EVERY SIX MONTHS		Check that there is no leaking of hydraulic oil or wear and tear of moveable components behind the lower cover of the cassette. Remove the lower cover to perform the inspection.
	I	Check that the two hinges of the hand rails are ***lubricated (I).
	J	Check that the four hinges, which connect the foldable parts of the platform, are ***lubricated (J).
	K	Check that the three bushings, holding the linkage arm (K) to the bridge plate, are ***lubricated.
	L	Check that all four support rollers are not deformed or damaged.
	M	Clean and **lubricate the tracks (M).
	N	Clean and ***lubricate the tracks (N). Reinstall the lower cover.
	O	Check that the two hinges (O) of the front hatch are ***lubricated and that the hatch is closing properly.
		Check that all bracket bolts are properly tightened.
	Check that the emergency release functions.	
	Check that bolts, pins and joints for wear/damage and that they are properly secured.	


There is an oil filter on the bottom of the oil container

## Technical specification

### Product data

Cassette length.....	1220 mm
Cassette width.....	1050 mm
Cassette height.....	140 mm
Usable platform length.....	1210-1340 mm
Usable platform width.....	800 mm
Minimum overlap, vehicle floor.....	20 mm
Minimum width, door opening.....	820 mm
Cassette weight.....	130 kg
Noise level during use.....	less than 70 dB

Pump module length.....	245 mm
Pump module width.....	120 mm
Pump module height.....	453 mm

### Product performance

Moving, vertical cycle.....	1075 mm
Moving from upper edge of cassette and up.....	430 mm
Minimum move up.....	130 mm
Moving from upper edge of cassette and down.....	645 mm
Load capacity.....	350 kg

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






- Connector/Terminal
- Diode (1.5KE18CA)
- Hydraulic Pump
- Relay
- Hydraulic Valve
- Motor Lift In/Out
- Relay
- Limit Switch
- Connection
- Signal Out (+12V) Lift Home Position
- Home Position +12V Lamp



Hydraulic diagram

Exploded view, pump module


30	1	Hydraulpump	438727	Vibo Power Unit CMB3-VB Motor: MD12080AN
29	2	Velcro 20mm, soft black	452172	452172
28	2	Velcro 20mm, tag black	452171	452171
27	2	Isolerande, Flexibelt lock SR 1044	400895	FCZF-0138-013
26	1	Gångad plugg med O-ring BSP 3/8x19	400875	400875
25	1	Kabelsko Stift 6,3x0,8 Röd	400859	Red 6,3x0,8 Tab 140896-2
24	1	Ring Terminal, PIDC, Red M5	411551	Red Ring M5 2-36153-4
23	2	Hex bolt 5/16-18 UNC - 0,75	440621	440621
22	1	Gummklämma	413141	G Klämna D8
21	1	Slang PA-10/8-100m N Längd 270	440450	Slang Ø10x1
20	1	Push-in rak koppling	440451	PLGA1006
19	1	Push-in rak koppl G3/8"O-ring 10 mm	440602	440602
18	1	Ändskydd	440614	15005
17	1	Ground Cable	440385	440385
16	1	Ground Cable	440384	440384
15	1	Power Peak Diode	440236	440236
14	2	Screw M8 x 12 ISO 4017 FZB 8.8	413143	413143
13	2	Serrated Lock Washer DIN 6798 A 8,4	440623	440623
12	1	Hand Control	440605	M102-3718
11	1	MÄTNIPPEL Ø8L MED MUTTER	440463	MD1-Ø8L-DKOL
10	1	T-adaptler 8L T Skarv	440462	XT08L
9	1	Adapter rak INV 8L - UTV. 1/4" NPT	440461	EGEØ8LN
8	1	Käpa Pumpfus	440133	M102-3709
7	1	Kontaktus plåt	424844	M102-3707
6	7	Relä	410806	NVF4-1C-S30AR
5	8	Screw M5 x 8 ISO 4762 FZBk 10.9	423464	423464
4	4	Plain washer 5,3x10x1 ISO 7089	413081	413081
3	1	Stag pumphus konsol	424843	M102-3702
2	1	Wire Harness Pump	452730	452730
1	1	Powerpack bracket	452005	452005
Det-nr.	Antal	Benämning	Art.nr.	Anmärkning
Design	Drawn	Checked	Project	Date
SER	ser	ATE	M102	2018-01-22
Power pack unit, Assy			Part no.	Scale
Mec			452631	1:2
Cassette lift			Revision	A

Autoadapt

Exploded view, Cassette lift

40	1	Adapter w/Trans Diode	451121	451121
39	1	Retaining Washer 8 DIN 6799 X.X	418158	418158
38	2	Washer, polyamid 10,2x20x1	452639	452639
37	4	Screw M5 x 8 ISO 7380 SS A2	440627	440627
36	1	Cross bar	452459	452459
35	1	Lever	452412	452412
34	1	Lever, Weld Assy	452621	452621
33	1	Outer platform, Black Assy	452622	452622
32	1	Screw M5 x 12 ISO 7380 FZBk 10.9	417756	417756
31	1	Packaging & fictitious parts	452875	452875
30	2	Tube ends	440617	440617
29	1	Ball bearing	440629	lbec 66H 625zz
28	1	Power pack unit, Assy	452631	452631
27	1	Cable sleeving Twist-in 25 L=1 M	440615	Kabelstrumpa 1M
26	1	Wire Harness in/Out Motor	452725	452725
25	1	Wire Loop On/Off	438792	M102-3715
24	1	Wire Harness Lift Switch	438791	M102-3710
23	1	Wire Harness Main Lift	452729	452729
22	2	Hose length 390 3/16	440448	Hydraulsiang Auto3
21	1	Hose length 1100 3/16	440447	Hydraulsiang Auto2
20	1	Hose length 2500 1/4	440446	Hydraulsiang Auto1
19	2	Angle adapter JIC 7/16-UNF 7/16 0.8	440458	JW/EEU-04-04-IDO,8
18	1	Hex Nut M8x1,25 ISO 4032 FZBk 8.8	418348	418348
17	5	Retaining Washer 5 DIN 6799 A2	438909	438909
16	1	Gas spring	438069	Gas-Spring-15-6EW-9075
15	1	Pivot pin	438725	M102-4060
14	1	Link	424966	M102-4042
13	1	Outer rod	438724	M102-4041
12	1	Inner rod	438723	M102-4040
11	2	Handle	424910	M102-4709
10	2	Screw Set M6x16 ISO 4029 SS A2	438722	438722
9	2	Hydraulic lift cylinder	440633	XE619300000HM83
8	1	Hatch, Mounting Assy	400828	M102-9108
7	1	Rollstop, Mounting Assy	452625	452625
6	1	Bridge plate, Black Assy	452087	452087
5	1	Inner platform, Black Assy	452085	452085
4	2	Parallel arm rod, Mounting Assy	400833	M102-9103
3	1	Lift arms, Mounting Assy	400834	M102-9102
2	1	Carriage, Mounting Assy	452624	452624
1	1	Cassette, Mounting Assy	400836	M102-9100
Det-nr	Antal	Benämning	Art.nr.	Anmärkning
Design	Drawn	Checked	Project	Date
Ser	Ser		M102	2018-03-22
Lift unit, Top Assy			Part no.	Scale
L1340 x W825 Mec Black ed.			452629	1:1
Cassette lift				Revision
				A

## Product certificate

### Manufacturer /

BraunAbility Europe AB  
Åkerivägen 7  
443 61 Stenkullen  
Sweden

Phone: +46 302 254 00  
E-mail: [info@braunability.eu](mailto:info@braunability.eu)  
[www.braunability.eu](http://www.braunability.eu)

### Declares that the product /

BraunAbility Cassette lift

### Conforms to following directives, standards and regulations /

2006/42/EEC  
EN 1756-2:2004+A1 2009  
UN ECE R-10 Rev 4:2012  
including amendment 2:2013

BraunAbility Europe AB, Stenkullen  
2019-04-02

Managing Director

Accredited by SWEDAC as testing laboratory in accordance with ISO/IEC 17025:2005  
Quality system certified in accordance with ISO 9001:2015

Illustrations, descriptions and specifications in the user manual are based on current product information.  
BraunAbility Europe AB reserves the right to make alterations without previous notice.  
© 2019 BraunAbility Europe AB