

# Occupant Restraints

User Manual



EN Instructions for fitting and use


Thank you for choosing

# Occupant restraints from BraunAbility!

The following manual is an important part of the product, providing you with information on how to achieve maximum performance and safe operation. Keep the manual in a safe place so that you can refer to it when necessary.

If you have any questions about your equipment, please contact us.

Once again, thank you for placing your confidence in our products!



Safe vehicle adaptation solutions  
*For your safety BraunAbility products are  
designed and tested according to current  
directives and standards.*

## Safety information

---

- 4 Limitations of use
- 6 General guidance
- 8 Before fitting and using the occupant restraints

## Fit and use

---

- 10 Fitting & using 3 point retractable occupant restraint Fitting
- 14 Fitting & using the 3 point static occupant restraint

## After care

---

- 20 Equipment storage and maintenance
- 20 Warranty

“Declaration of conformity” at the end of the manual.

# Safety information

## Limitation of use

The development of BraunAbility wheelchair and occupant restraints is a continuous process. Applications are added on a regular basis. For more details contact BraunAbility or look on the BraunAbility website.



*These products have been designed to secure the passenger only and must only be used where the wheelchair has already been secured with a separate wheelchair tie-down system.*



*In addition to the details given in these instructions, users of BraunAbility restraints must refer to the wheelchair manufacturer's 'Instructions for Use in Transport' for full details of tie-down attachment points on the wheelchair, plus any other specific instructions relating to use in transport.*



*Wheelchair users and their carers must make sure that their wheelchair is recommended for use in transport, including any 'add-on' components such as power tilt or recline options.*



*BraunAbility recommends the use of a 3 point occupant restraint system to provide greater protection in case of an impact. However we also recognise that some vehicle layout/designs or specific postural or medical conditions do not allow that style of seatbelts to be used easily.*

*This system is not ISO 10542 compliant if used in 2 point configuration.*



*In accordance with international regulations, our occupant restraints are designed to be used for passengers weighing 22kg or above. If the passenger weighs less than 22kg, we recommend that a suitable, and appropriately tested child restraint seat is used. This may involve a secondary seat belt restraint as recommended by the original seat manufacturer.*

## Occupant Restraints

OCR01 - 3 Point Retractable  
Occupant Restraint



OCR02 - 3 Point Static  
Occupant Restraint

## General guidance

- Wheelchair Accessories that have not been approved by the Wheelchair Manufacturer must be removed from the wheelchair and secured in the vehicle during transport to reduce the potential for injury. Refer to 'Instructions For Use in Transport' provided with the wheelchair or contact wheelchair manufacturer for further guidance.
- These occupant restraints comply with all applicable requirements of ISO 10542, including a 48km/h, 20g frontal impact test using a forward facing 85kg surrogate wheelchair and an ATD (test dummy) with a mass of 76.3kg. The test dummy was restrained by both a pelvic and upper-torso restraint. Use of a pelvic only belt may compromise the performance of the WTORS system and should be avoided.
- If the installation is to be used with an occupant head rest anchored to the vehicle, then a vehicle-anchored back rest must be provided to minimise rearward deflection of the wheelchair seatback, preventing neck injury.
- Regular inspection of all parts is recommended and the equipment should be used only if all components are in good condition.
- **Warning:** protect webbing from contacting sharp edges and corners. Replace equipment if the webbing does becomes cut, contaminated or frayed.

- Any restraints which have been subjected to a crash situation from which the vehicle must be towed, should, in the interest of safety, be replaced.
- Do not attempt to modify the equipment. For further advice on the installation and use of this tie-down, please contact BraunAbility. Please read this manual fully before use.
- Avoid contact with corrosive liquids. Care should be taken to prevent contamination of the webbing with polishes, oils and chemicals.
- The tie-down anchorages should be installed by an experienced technician/vehicle converter. Anchorages should not be installed into unsound materials such as corroded metal, wood, plastic and fibre glass panels, without additional and suitable reinforcement.
- The equipment has been tested in a configuration similar to that recommended by BraunAbility, and any deviation from the recommendations here is the responsibility of the installer/user.

## Before fitting and using the occupant restraints

- Ensure that the wheelchair is correctly maintained and that the settings of any adjustable parts are made according to 'Instructions for Use in Transport'.
- Whenever possible remove any items of luggage etc that may be attached to the wheelchair and secure or store separately during transport in order to reduce the potential for injury to other passengers travelling in the vehicle.
- Extra care must be exercised when using vehicle boarding aids such as passenger lifts or ramps during the loading process. Refer to 'Instructions For Use' for information on safe slopes.
- Position the wheelchair facing forward - centrally in the designated region of the vehicle.  
**Ensure the wheelchair brakes are applied.**
- To minimise the potential for head injuries in an impact, allow a clear space of at least 400 mm behind and 650 mm, (FCZ, front clear zone), in front of the head of the wheelchair user, (Fig A). The shoulder belt anchorage must be roof or side-wall anchored at a height level such that the belt webbing passes over the midpoint of the occupant shoulder and at a height that is at or just above the level of the occupants shoulders so as not to impose downward loads on the spine.
- A height provision (HHT) ranging from 1000mm to 1550mm should be made, depending on the size of the passenger. There should also be 200 mm of clear space either side of the wheelchair centre line. If these clear space dimensions cannot be provided then any should be adequately padded and comply with impact performance requirements of ECE Regulation 21 'FMVSS 201'. All vehicle padding should comply with the flammability requirements of ECE Regulation 118 'FMVSS 302'.

*Note: seated head height (HHT) ranges from as low as 1000mm for a 6-year-old child to 1550mm for a tall adult.*



- Wheelchair users, their carers and family are advised to check vehicle specifications to ensure that sufficient floor space is available to accommodate the wheelchair and tie-down system. These distances are based upon the desire to maintain clear zones for potential head excursions of occupants provided with both upper and lower torso restraints.
- Users of heavy powered wheelchairs are also advised to check vehicle carrying capacity. If in doubt consult the vehicle supplier for further details.
- Any airbag, as fitted to the vehicle, shall be used only as a supplementary occupant restraint if designed to be used in combination with the wheelchair tie-down.
- Installers of this tie-down should take note of any vehicle airbag position when planning the installation. Airbags can cause serious injury if a wheelchair-seated occupant is seated too close to an airbag position. If in doubt contact the vehicle manufacturer or your National Automotive Regulatory Body for advice.




Figure A

# Fit and use

## Fitting and using the 3 point retractable occupant restraint

### Rail Floor

1. The rail will have been installed in the vehicle in accordance with our own and the vehicle converter's instruction. Position wheelchair within vehicle as required.
2. Wheelchair tie-downs (not supplied with these restraints) should be fitted first to secure the wheelchair, before any occupant restraint is fitted.
3. Position the occupant restraint behind the wheelchair tie-down, with reel housing (Fig. 1A) on the window side of the vehicle and reversible stalk (Fig 1B) on the aisle side.
4. Fit the occupant restraint ATF (aluminium track fitting) into the rail by aligning the ATF feet with the cut-out sections of the rail. Locate into the rail, (Fig 2).
5. Press on the ribbed part of the ATF, (Fig 2A), and push firmly towards the wheelchair until the yellow plunger drops and locks into the rail.
6. Position the third point fixing into the cant rail (Fig 3) so that it is vertically above the inertia reel casing fixed onto the floor fixing. Unfasten the tongue and buckle.
7. Pull the webbing upwards and unfasten the tongue from the red webbing buckle.
8. Position the black webbing to form the lap belt and insert the tongue into the buckle stalk. Ensure that the lap belt lays low on the pelvis of the occupant, running as close as possible over the hips on both sides.
9. Remove the black plastic cover on the red shoulder belt tongue and fit into the third point fixing. Position the red webbing to form the shoulder belt and insert the buckle into the tongue. Adjust the height of the shoulder belt to clear the occupant's shoulder by approximately 25mm/1 inch (Fig 4).
10. The lap belt anchor points should be positioned to achieve belt angles of 30° or more to the horizontal and preferably between 45° and 75° in order to fit low across the pelvis reducing the possibility of the belt loading the abdomen (Fig 5). The pelvic restraint is designed to bear upon the bony structure of the body and should be worn low across the front of the pelvis with any junctions between the pelvic and shoulder restraints located near the wearers hips.



### Removing the occupant restraint

1. Release the red shoulder belt from the lap belt section, then from the third point fixing and carefully allow the shoulder belt to return to the retractor.
2. Disconnect the black lap belt from the buckle stalk on the aisle side and carefully allow the lap belt webbing to return onto its retractor.
3. Connect the shoulder belt buckle to the lap belt tongue above the opening of the retractor box.
4. The wheelchair tie-down must now be removed.




Figure 2




Figure 3



Figure 4




Figure 5



Figure 1

# Fit and use

## Fitting and using the 2 point retractable occupant restraint

### Rail floor

1. The rail will have been installed in the vehicle in accordance with our own and the vehicle converter's instruction. Position wheelchair within vehicle as required.
2. Wheelchair tie-downs (not supplied with these restraints) should be fitted first to secure the wheelchair, before any occupant restraint is fitted.
3. Position the occupant restraint behind the wheelchair tie-down, with reel housing (Fig. 1A) on the window side of the vehicle and reversible stalk (Fig 1B) on the aisle side.
4. Fit the occupant restraint ATF (aluminium track fitting) into the rail by aligning the ATF feet with the cut-out sections of the rail. Locate into the rail, (Fig 2).
5. Press on the ribbed part of the ATF, (Fig 2A), and push firmly towards the wheelchair until the yellow plunger drops and locks into the rail.
6. Ensure the Comfort Clip is released (Fig 7) and press the buckle to release the lap belt.
7. Pull the metal tongue on the black webbing upwards to release the webbing from both inertia reels. Position the webbing across the occupant with the black webbing forming the lap belt and the red webbing forming the shoulder belt (Fig 6).
8. Insert the tongue fixed to the black webbing into the stalk buckle.
9. Adjust for comfort by raising the webbing on the shoulder by two finger widths, unlock the Comfort Clip to take up the slack, then re-lock the clip so that it rests on the reel cover.
10. The lap belt anchor points should be positioned to achieve belt angles of 30° or more to the horizontal and preferably between 45° and 75° in order to fit low across the pelvis reducing the possibility of the belt loading the abdomen (Fig 5). The pelvic restraint is designed to bear upon the bony structure of the body and should be worn low across the front of the pelvis with any junctions between the pelvic and shoulder restraints located near the wearers hips.



### Removing the occupant restraint

1. Release the Comfort Clip.
2. Release the red shoulder belt from the lap belt section, and carefully allow the shoulder belt to return to the retractor.
3. Disconnect the black lap belt from the buckle stalk on the aisle side and carefully allow the lap belt webbing to return onto its retractor.
4. Connect the shoulder belt buckle to the lap belt tongue above the opening of the retractor box.
5. The wheelchair tie-down must now be removed.




Figure 5




Figure 2



Figure 7




Figure 1



Figure 6

# Fit and use

## Fitting and using the 3 point static occupant restraint

### Rail floor

1. The rail will have been installed in the vehicle in accordance with our own and the vehicle converter's instruction. Position wheelchair within vehicle as required.
2. Wheelchair tie-downs should be fitted first to secure the wheelchair, before any occupant restraints are fitted.
3. Position the black webbing lap belt (A1 & A2) with the ATF feet facing the towards the back of the wheelchair and the webbing side with the male pin connector (A2) located on the opposite side of the third point fixing wall.
4. Fit the lap belt track fitting into the rail by aligning the feet with the cut-out sections of the rail and locate into the rail (Fig 2). Press on the ribbed part of the track fitting, (Fig 2), and push firmly towards the wheelchair until the yellow plunger drops and locks into the rail. Repeat operation with other side. **IMPORTANT:** Check plungers are fully engaged on both sides.
5. Position the lap belt across the occupant waist, clip the lap belt tongue and buckle together to fasten the lap belt around the occupant waist.
6. Position the third point fixing (C) into the cant rail (Fig 3), unfasten the tongue and buckle.
7. Fit the red webbing shoulder belt (B) track fitting into the rail by aligning the feet with the cut-out sections of the rail and locate into the rail. Press on the ribbed part of the track fitting (Fig 2) and push firmly towards the wheelchair until the yellow plunger drops and locks into the rail. **IMPORTANT:** Check plunger is fully engaged.
8. Position the red webbing to form the shoulder belt and insert the buckle into the 3rd point tongue (Fig 9).
9. Fasten the female pin connector on the shoulder belt (Fig 10) to the male pin connector on the waist belt to form the shoulder belt.
10. Adjust the height of the shoulder belt to clear the occupant's shoulder by approximately 25mm (1 inch) Fig 4.
11. The lap belt anchor points should be positioned to achieve belt angles of 30° or more to the horizontal and preferably between 45° and 75° in order to fit low across the pelvis reducing the possibility of the belt loading the abdomen (Fig 5). The pelvic restraint is designed to bear upon the bony structure of the body and should be worn low across the front of the pelvis with any junctions between the pelvic and shoulder restraints located near the wearers hips.



### Removing the occupant restraint

1. Unfasten the third point from the shoulder belt and reconnect the tongue and buckle on the third point fixing.
2. Lift the yellow retaining rail fixing clip fully and remove the third point fixing from the wall rail.
3. Disconnect the pin connector and remove the shoulder belt from the passenger's shoulder.
4. Lift the yellow retaining rail fixing clip fully and remove the shoulder belt from the rail.
5. Unfasten the lap belt tongue and buckle and remove the lap belt from the occupant.
6. Lift the yellow retaining rail fixing clips fully and store the occupant restraint safely.



Figure 8




Figure 2



Figure 3

EN

Fit and use



Figure 9



Figure 10



Figure 4




Figure 5

# After care

## Equipment storage and maintenance

- Store the restraint safely off the floor to avoid damage and ensure that it cannot become a projectile in an accident. This can be achieved by using wall mounted storage bag such as SLR111.
- Regularly inspect the restraint systems for damage, wear or malfunction. If any problems are identified replace it immediately.
- When not in use, keep loose occupant webbing ends connected to their corresponding buckle sleeves, etc, to prevent them from becoming tripping hazards and from flailing around when the vehicle is in motion.
- All webbing and components can be cleaned as necessary, but care should be taken to prevent contamination of the webbings with polishes, oils and chemicals, particularly battery acid.
- To clean the straps use warm soapy water and a clean soft cloth. Rinse with clear water and allow to air dry. To disinfect, use a mild spray disinfectant and do not use products containing bleach. **Important: when cleaning or disinfecting, do not immerse or flood buckles, karabiners fittings or floor anchors in the disinfectant or water.**
- If the vehicle is involved in an accident when any restraints are deployed, remove them from service and replace immediately. If in doubt please contact BraunAbility.



## After care

### Warranty

BraunAbility products are extensively tested using BraunAbility anchorage systems, and our full warranty normally only applies to BraunAbility equipment when used with BraunAbility branded anchorages or as instructed. BraunAbility have also participated in test programs with other manufacturers anchorage products and will support warranty on the BraunAbility products when used in conjunction with such jointly tested systems. For further details on specific applications please contact the Sales Office. In other situations, using BraunAbility products, for which BraunAbility has not participated in a joint test program, a limited BraunAbility warranty will apply.

# Declaration of conformity

## **Manufacturer /**

BraunAbility UK Ltd.  
Unwin House  
The Horseshoe Coat Road  
Martock, Somerset,  
TA12 6EY, UK

Phone: +44 1935 827740  
E-mail: sales.uk@braunability.eu  
www.braunability.eu

## **Declares that the products /**

OCR01, OCR02

## **Conforms to following directives standards and regulations /**

214/2014/EU Paragraphs 2.3.1 and 2.3.2  
ISO 10542:2012  
RESNA WC-4 2012 Section 18  
R.107.06  
2007/46/EC

BraunAbility, Martock, 14 March 2019

A handwritten signature in black ink, appearing to read 'Rob Butcher', with a long horizontal stroke extending to the right.

Rob Butcher  
Director of Engineering



Accredited by URS as testing laboratory in accordance with ISO/IEC 17025:2005  
Quality system certified in accordance with ISO 9001:2008

Illustrations, descriptions and specifications in the user manual are based on current product information. BraunAbility UK Ltd reserves the right to make alterations without previous notice.  
© 2019 BraunAbility UK Ltd

**BraunAbility UK Ltd**  
Unwin House  
The Horseshoe, Coat Road  
Martock, Somerset, TA12 6EY  
UK

Phone: +44 1935 827 740  
E-mail: [info@braunability.eu](mailto:info@braunability.eu)  
[www.braunability.eu](http://www.braunability.eu)